



F8 Antibody

Product Code	CSB-PA007932GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00451
Immunogen	Human F8
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IF
Form	Liquid
Storage Buffer	PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	F8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.