



VASH1 Antibody

Product Code	CSB-PA745337GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L8A9
Immunogen	Human VASH1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Form	Liquid
Storage Buffer	PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	VASH1 KIAA1036 VASH
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Cycle
Target Names	VASH1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.