




Recombinant Human Nicotinamide mononucleotide adenylyltransferase 3 (NMNAT3)

Product Code	CSB-EP839418HU
Relevance	Catalyzes the formation of NAD ⁺ from nicotinamide mononucleotide (NMN) and ATP. Can also use the deamidated form; nicotinic acid mononucleotide (NaMN) as substrate with the same efficiency. Can use triazofurin monophosphate (TrMP) as substrate. Can also use GTP and ITP as nucleotide donors. Also catalyzes the reverse reaction, i.e. the pyrophosphorolytic cleavage of NAD ⁺ . For the pyrophosphorolytic activity, can use NAD ⁺ , NADH, NaAD, nicotinic acid adenine dinucleotide phosphate (NAD ⁺), nicotinamide guanine dinucleotide (NGD) as substrates. Fails to cleave phosphorylated dinucleotides NADP ⁺ , NADPH and NaADP ⁺ . Protects against axonal degeneration following injury.
Abbreviation	Recombinant Human NMNAT3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96T66
Alias	Nicotinamide-nucleotide adenylyltransferase 3 ;NMN adenylyltransferase 3;Nicotinate-nucleotide adenylyltransferase 3 (EC:2.7.7.18) ;NaMN adenylyltransferase 3;Pyridine nucleotide adenylyltransferase 3 (EC:2.7.7.1) ;PNAT-3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MYQVIQGIISPVNDTYGKKDLAASHHRVAMARLALQTSDWIRVDPWESEQAQ WMETVKVLRHHHSHKLLRSPQMEGPDHGKALFSTPAAVPELKLLCGADVLT FQTPNLWKDAHIQEIVEKFGLVCVGRVGHDPKGYIAESPILRMHQHNIHLAKEP VQNEISATYIRRALGQGQSVKYLIPDAVITYIKDHGLYTKGSTWKGKSTQSTEG KTS
Research Area	Metabolism
Source	E.coli
Target Names	NMNAT3
Expression Region	1-215aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	40.1kDa

**Protein Length**

Full Length of Isoform 2

Image

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.