




Recombinant Human Tektin-5 (TEKT5)

Product Code	CSB-EP856979HU
Relevance	Structural component of ciliary and flagellar microtubules. Forms filamentous polymers in the walls of ciliary and flagellar microtubules .
Abbreviation	Recombinant Human TEKT5 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96M29
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MEFLGTTQTASYCGPKKCCGLTSLPAVQAPVIQECYQPYYPGYRYLNSWRP SLFYKIANVQTCPESTSTLRPPTILPTLRSALFSRYSPhdWDQSNQLQVRGA EASRLWASRLTDDSMRLQLQDKDQLTHQMQEGTCRNLGQRLSDIGFWKSELS YELDRLLTENQNLETVKRRLECAANEVNCPLQVALECLYHREKRIGIDLVDHNV EKNLIREVDLLKCCQEQRKLAQRIDIQMRDNRDAQHVLERDLEDKSSAQCID EKCFNLRNTSDCISFFHGMKIDGTISVPETWAKFSNDNIKHSQNMRRANSIQLR EEAEHLFETLSDQMWRQFTDTNLAFNARISEVTDVKNKLQTQLAKTLQEIFA ENTIMLLERSIMAKEGPKVAQTRLECRTRRPNMELCRDIPQLKLVNEVFTIDD TLQTLKLRLRETQDTLQLLVMTKCRLEHELAIKANTLCIDKEKCMGMRKTFPCT PRLVGHT
Research Area	Signal Transduction
Source	E.coli
Target Names	TEKT5
Expression Region	1-485aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	72.3kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.