




# Recombinant Human N-acetylneuraminase 9-O-acetyltransferase (CASD1), partial

<b>Product Code</b>	CSB-EP846683HU
<b>Abbreviation</b>	Recombinant Human CASD1 protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96PB1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	ASRRYRGNDSC EYLLSSGRFLG EKVWQPHSCMMHKYKISEAKNCLVDKHIAFI GDSRIRQLFY SFVKIINPQFKEEGNKHENIPFEDKTASVKVDFLWHPEVNGSMK QCIKVVWTE DSIKPHVIVAGAATWSIKIHNGSSEALSQYKMNITSIAPLLEKLAKT SDVYWVLQDPVYEDLLSEN RKMITNEKIDAYNEAAVSILNSSTRNSKSNVKMF SVSKLIAQETIMESLDGLHLP ESSRETTAMILMNVYCNKILKPVDGSCCQPRPP VTLIQ
<b>Research Area</b>	Cell Biology
<b>Source</b>	E.coli
<b>Target Names</b>	CASD1
<b>Expression Region</b>	39-313aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-SUMO-tagged
<b>Mol. Weight</b>	47.2kDa
<b>Protein Length</b>	Partial

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.