




Recombinant Human Ribonuclease P protein subunit p29 (POP4)

Product Code	CSB-EP018373HU
Relevance	Part of ribonuclease P, a protein complex that generates mature tRNA molecules by cleaving their 5'-ends. May function with RPP38 to coordinate the nucleolar targeting and/or assembly of RNase P.
Abbreviation	Recombinant Human POP4 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O95707
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MKSVIYHALSQKEANDSDVQPSGAQRAEAFVRAFLKRSTPRMSPQAREDQLQ RKAVVLEYFTRHKRKEKKKKAKGLSARQRRELRLFDIKPEQQRYSLFLPLHEL WKQYIRDLC SGLKPDTQPQMIQAKLLKADLHGAIISVTKSKCPSYVGITGILLQE TKHIFKIITKEDRLKVIPKLNCFVTVETDGFISYIYGSKFQLRSSERSAKKFKAKG TIDL
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	POP4
Expression Region	1-220aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	41.4kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.