



Recombinant Human DNA repair protein XRCC3 (XRCC3)


Product Code	CSB-EP026231HU
Relevance	Involved in the homologous recombination repair (HRR) pathway of double-stranded DNA, thought to repair chromosomal fragmentation, translocations and deletions. Part of the RAD21 paralog protein complex CX3 which acts in the BRCA1-BRCA2-dependent HR pathway. Upon DNA damage, CX3 acts downstream of RAD51 recruitment; the complex binds predominantly to the intersection of the four duplex arms of the Holliday junction (HJ) and to junctions of replication forks. Involved in HJ resolution and thus in processing HR intermediates late in the DNA repair process; the function may be linked to the CX3 complex and ses to involve GEN1 during mitotic cell cycle progression. Part of a PALB2-scaffolded HR complex containing BRCA2 and RAD51C and which is thought to play a role in DNA repair by HR. Plays a role in regulating mitochondrial DNA copy number under conditions of oxidative stress in the presence of RAD51 and RAD51C.
Abbreviation	Recombinant Human XRCC3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O43542
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MDLDLLDLNPRIIAAIKKAKLKSVEVLHFSGPDLKRLTNLSSPEVWHLLRTASL HLRGSSILTALQLHQKQERFPTQHQRSLGCPVLDALLRGGLPLDGITELAGRS SAGKTQLALQLCLAVQFPRQHGGLGAGAVYICTEDAFPHKRLQQLMAQQPRL RTDVPGELLQKLRFSGSQIFIEHVADVDTLLECVNKKVPVLLSRGMARLVVIDSV AAPFRCEFDSQASAPRARHLQSLGATLRELSSAFQSPVLCINQVTEAMEEQGA AHGPLGFWDERVSPALGITWANQLLVRLADRLREEEAALGCPARTLRVLSAP HLPPSSCSYTISAEGVRGTPGTQSH
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	XRCC3
Expression Region	1-346aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged



Mol. Weight 50.8kDa

Protein Length Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.