




Recombinant Human Coiled-coil domain-containing protein 36 (CCDC36), partial

Product Code	CSB-EP816910HU
Abbreviation	Recombinant Human CCDC36 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IYA8
Alias	Cancer/testis antigen 74 ;CT74
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MLSIPSGSGNKKSSNWNNNQNDYSSLSDSQFLFGSQFCPENSETLSAPLDFG AHLRHSKQSQQNYLEGEPSIFTKYQTKPQLFGGDIKDGGFLPPPLSVGKSKGL LEQFEEKKKRAKDKCDSETLYNFVSNVRESILRLQTSVEKSEDHLSSRSQSILD SLETVAKTLQETIQAQNDLVFEAVQDKGNMEQAILEMKKRFEARQGEFIEMKS NLKHLEVLVAQQSQEFQQLCEQLGQLNVPSVLAELKRLISVPPVKDSASQTSP PLAQSLNLTRQEKYTSEKPVLWQAQALPAAWNPGMGS LQPGEFDVWGEGAK NDDLQEEAALPAFGSHERNRHVKDKVVQTNCKNWAVTKTGAKNHGSSVPGH KIPSDRDLVSQGASQLTSLEINFSTSIKNACQKYQAQSMFLCDPREHLVIKQKD GTVEMRGKDKKQQPRKAHRAHRGRLIASKQKQIPIQTCKFNSKYQSPQPAISV PQSPFLGQQEPRAQPLHLQCPRSPRPVCPILGGTVMPNKT VRAVQGRLLQL SRCSSQDNWLLSSSSQGDHQMWSFSDLNLGCSETPLCKEAGKNLLYDLGFD SSDDD
Research Area	Cell Biology
Source	E.coli
Target Names	CCDC36
Expression Region	11-592aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	80.9kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.