




Recombinant Escherichia coli Long-chain-fatty-acid--CoA ligase (fadD)

Product Code	CSB-EP300208ENV
Relevance	Catalyzes the esterification, concomitant with transport, of exogenous long-chain fatty acids into metabolically active CoA thioesters for subsequent degradation or incorporation into phospholipids. Activity is the highest with fatty acid substrates of > 10 carbon atoms
Abbreviation	Recombinant E.coli fadD protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P69451
Alias	Long-chain acyl-CoA synthetase Short name:Acyl-CoA synthetase
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MKKVWLNRYPADVPTEINPDRYQSLVDMFEQSVARYADQPAFVNMGEVMTF RKLEERSRAFAAYLQQGLGLKKGDRVALMMPNLLQYPVALFGILRAGMIVVNV NPLYTPRELEHQLNDSGASAIVIVSNFAHTLEKVVDKTAVQHVILTRMGDQLST AKGTVVNFVVKYIKRLVPKYHLPDAISFRSALHNGYRMQYVKPELVPEDLAFLO YTGGTTGVAKGAMLTHRNMLANLEQVNATYGPLLHPGKELVVTALPLYHIFAL TINCLLFIELGGQNLLITNPRDIPGLVKELAKYPFTAITGVNTLNFALLNNKEFQQ LDFSSLHLSAGGGMPVQQVVAERWVKLTGQYLLEGYGLTECAPLVSVPYDI DYHSGSIGLPVPSTEAKLVDDDDNEVPPGQPGELCVKGPQVMLGYWQRPDA TDEIIKNGWLHTGDIAVMDEEGFLRIVDRKKDMILVSGFNVYPNEIEDVVMQHP GVQEVAAVGVPSGSSGEAVKIFVVKKDPSSLTEESLVTFCRRQLTGYKVPKLVE FRDELPKSNVGKILRREL RDEARGKVDNKA
Research Area	Others
Source	E.coli
Target Names	fadD
Protein Names	Long-chain acyl-CoA synthetase Short name:Acyl-CoA synthetase
Expression Region	1-561aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	66.3 kDa

**Protein Length**

Full Length

Image

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.