




Recombinant Mouse Apolipoprotein E (ApoE)

Product Code	CSB-EP001936MO
Relevance	Mediates the binding, internalization, and catabolism of lipoprotein particles. It can serve as a ligand for the LDL (apo B/E) receptor and for the specific apo-E receptor (chylomicron remnant) of hepatic tissues.
Abbreviation	Recombinant Mouse ApoE protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P08226
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	EGEPEVTDQLEWQSNQPWEQALNRFWDYLRWVQTLSDQVQEELQSSQVTQ ELTALMEDTMTEVKAYKKELEEQLGPVAEETRARLGKEVQAAQARLGADMED LRNRLGQYRNEVHTMLGQSTEEIRARLSTHLRKMRLMRDAEDLQKRLAVY KAGAREGAERGVS AIRERLGPLVEQGRQRTANLGAGAAQPLRDRAQAFGDRI RGRLEEVGNQARDRLEEVRHEMEEVRSKMEEQTQQIRLQAEIFQARLKGWFE PIVEDMHRQWANLMEKIQASVATNPITPVAQENQ
Research Area	Others
Source	E.coli
Target Names	ApoE
Expression Region	19-311aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	38.0kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.