



Recombinant Mouse Complement C3 (C3), partial

Product Code	CSB-MP365507MO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P01027
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	SEETKQNEAFSLTAKGKGRGTLSSVVAVYHAKLKSQVTCKKFDLRVSIRPAPET AKKPEEAKNTMFLEICTKYLGDVDATMSILDISMGTGFAPDTKDLELLASGVDR YISKYEMNKAFSNKNTLIIYLEKISHTTEEDCLTFKVHQYFNVGLIQPGSVKVSYSY YNLEESCTRFYHPEKDDGMLSKLCHSEMCRCAEENCFMQQSQEKINLNVRLD KACEPGVDYVYKTELNIELLDDFDEYTMTIQQVIKSGSDEVQAGQQQRKFISHI KCRNALKLQKGKKYLMWGLSSDLWGEKPNTSYIIGKDTWVEHWPEAEECQD QKYQKQCEELGAFTESMVVYGCPN
Source	Mammalian cell
Target Names	C3
Protein Names	Recommended name: Complement C3 Alternative name(s): HSE-MSFC Cleaved into the following 10 chains: 1. Complement C3 beta chain 2. Complement C3 alpha chain 3. C3a anaphylatoxin 4. Complement C3b alpha chain 5. Complement C3c a
Expression Region	1321-1663aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.