




Recombinant Human B-lymphocyte antigen CD20 (MS4A1), partial

Product Code	CSB-YP015007HU1
Relevance	This protein may be involved in the regulation of B-cell activation and proliferation.
Abbreviation	Recombinant Human MS4A1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P11836
Alias	B-lymphocyte surface antigen B1Bp35Leukocyte surface antigen Leu-16Membrane-spanning 4-domains subfamily A member 1; CD20
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	IKISHFLKMESLNFIRAHTPYINIYNCEPANPSEKNSPSTQYCYSIQSLFLGILSV MLIFAFFQELVIAGIVENEWKRTC SRPKSNIVLLSAEEKKEQTIEIKEEVVGLTET SSQPKNEEDIEIPIQEEEEETETNFPEPPQDQESSPIENDSSP
Research Area	Immunology
Source	Yeast
Target Names	MS4A1
Expression Region	141-297aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	20.0kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.