




# Recombinant Bacillus subtilis Malate dehydrogenase (mdh)

<b>Product Code</b>	CSB-EP013623BRJ
<b>Relevance</b>	Catalyzes the reversible oxidation of malate to oxaloacetate.
<b>Abbreviation</b>	Recombinant Bacillus subtilis mdh protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P49814
<b>Alias</b>	Vegetative protein 69 ;VEG69
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bacillus subtilis (strain 168)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	GNTRKKVSVIGAGFTGATTAFLIAQKELADVVLVDIPQLENPTK GKALDMLEAS PVQGFDAKITGTSNYEDTAGSDIVVITAGIARKPGMSRDDLSTNEKIMRSVTQ EIVKYS PDSIIVL TNPVDAMTYAVYKESGFPKERVIGQSGVLD TARFRTFVAEE LNLSVKDVTG FVLGGHGDDMVPLV RYSYAGGI PLET LIPKERIDAIVERTRKGG GEIVNLLGNGSAYYAPAASLTEMVEAILKDQRRVLPTIAYLEGEYGYEGIYLGVP TIVGGNGLEQII ELELTDYERAQLNKSVESVKNVMKVLS
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	mdh
<b>Expression Region</b>	2-312aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	37.5kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.