




Recombinant Mouse H-2 class I histocompatibility antigen, K-D alpha chain (H2-K1), partial

Product Code	CSB-EP355780MO
Relevance	Involved in the presentation of foreign antigens to the immune syst.
Abbreviation	Recombinant Mouse H2-K1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P01902
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	GPHSLRYFVTAVSRPGLGEPRIAVGYVDDTQFVRFDSDADNPRFEPRAPWM EQEGPEYWEEQTQRAKSDEQWFRVSLRTAQRYYNQSKGGSHTFQRMFGCD VGSDWRLRLRGYQQFAYDGRDYIALNEDLKTWTAADTAALITRRKWEQAGDAE YYRAYLEGECEVWLRRLYLELGNETLLRTDSPKAHVITYHPRSQVDVTLRCWAL GFYPADITLTWQLNGEDLTQDMELVETRPAGDGTFFQKWAAVVPLGKEQNYT CHVHHKGLPEPLTLRWKLPPSTVSNT
Research Area	Others
Source	E.coli
Target Names	H2-K1
Expression Region	22-305aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	60.0kDa
Protein Length	Extracellular Domain
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.