




# Recombinant Human Protein-tyrosine sulfotransferase 2 (TPST2), partial

<b>Product Code</b>	CSB-EP024133HU
<b>Relevance</b>	Catalyzes the O-sulfation of tyrosine residues within acidic motifs of polypeptides.
<b>Abbreviation</b>	Recombinant Human TPST2 protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O60704
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	QQVLECRAVLAGLRSPRGAMRPEQEELVMVGTNHVEYRYGKAMPLIFVGGVP RSGTTLMRAMLDHAHPEVRCGEETRIIPRVLAMRQAWSKSGREKRLRDEAGVT DEVLDAAMQAFILEVIAKHGEPARVLCNKDPFTLKSSVYLSRLFPNSKFLLMVR DGRASVHSMITRKVTIAGFDLSSYRDCLTKWNKAIEVMYAQCMEVGKEKCLPV YYEQLVLHPRRSLKLILDFLGIAWSDAVLHHEDLIGKPGGVSLSKIERSTDQVIK PVNLEALSKWTGHIPGDVVRDMAQIAPMLAQLGYDPYANPPNYGNPDPFVINN TQRVLKGDYKTPANLKG YFQVNQNSTSSHLGSS
<b>Research Area</b>	Cardiovascular
<b>Source</b>	E.coli
<b>Target Names</b>	TPST2
<b>Expression Region</b>	26-377aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-SUMO-tagged
<b>Mol. Weight</b>	55.3kDa
<b>Protein Length</b>	partial
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.