



Recombinant Bovine Guanine nucleotide-binding protein subunit alpha-14 (GNA14)

Product Code	CSB-EP009586BO-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P38408
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	MAGCCCLSAE EKESQRISAE IERQLRRDKK DARRELKLLL LGTGESGKST FIKQMRIIHG SGYSDEDRKG FTKLVYQNI F TAMQAMIRAM DTLKIQYVCE QNKENAQLIR EVEVDKYSTL SRDQVEAIKQ LWQDPGIQEC YDRRREYQLS DSAKYYLTDI DRIAMPAFVP TQQDVLVRV PTTGIIIEYPF DLENIIFRMV DVGGQRSERR KWIHCFSVT SIIFLVALSE YDQVLAECDN ENRMEESKAL FKTIITYPWF LNSSVILFLN KKDILLEEKIM YSHLISYFPE YTGPKQDVKA ARDFILKLYQ DQNPDKKVI YSHFTCATDT ENIRFVFAAV KDTILQLNLR EFNLV
Source	E.coli
Target Names	GNA14
Protein Names	Recommended name: Guanine nucleotide-binding protein subunit alpha-14 Short name= G alpha-14 Short name= G-protein subunit alpha-14 Alternative name(s): G-protein subunit GL1 alpha
Expression Region	1-355
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.