




# Recombinant Human Neuropeptide FF receptor 2(NPFFR2)

|                          |                                                                                                                                                                                                                                                                                    |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-YP015983HU                                                                                                                                                                                                                                                                     |
| <b>Abbreviation</b>      | Recombinant Human NPFFR2 protein                                                                                                                                                                                                                                                   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| <b>Uniprot No.</b>       | A0PJM9                                                                                                                                                                                                                                                                             |
| <b>Product Type</b>      | Recombinant Proteins                                                                                                                                                                                                                                                               |
| <b>Immunogen Species</b> | Homo sapiens (Human)                                                                                                                                                                                                                                                               |
| <b>Purity</b>            | ≥ 85% as determined by SDS-PAGE.                                                                                                                                                                                                                                                   |
| <b>Research Area</b>     | Neuroscience                                                                                                                                                                                                                                                                       |
| <b>Source</b>            | Yeast                                                                                                                                                                                                                                                                              |
| <b>Target Names</b>      | NPFFR2                                                                                                                                                                                                                                                                             |
| <b>Tag Info</b>          | N-terminal His-tagged<br>Tag-Free                                                                                                                                                                                                                                                  |
| <b>Mol. Weight</b>       | 50.7 kDa                                                                                                                                                                                                                                                                           |
| <b>Protein Length</b>    | Full Length of Mature Protein                                                                                                                                                                                                                                                      |

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.  
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.