




Recombinant Mouse Serum amyloid A-3 protein (Saa3), partial

Product Code	CSB-EP361411MO
Relevance	Major acute phase reactant. Apolipoprotein of the HDL complex.
Abbreviation	Recombinant Mouse Saa3 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P04918
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	RWVQFMKEAGQGSRDMWRAYSDMKKANWKNSDKYFHARGNYDAARRGPG GAWAAKVISDAREAVQKFTGHGAEDSRADQFANEWGRSGKDPNHFRPAGLP KRY
Research Area	Others
Source	E.coli
Target Names	Saa3
Expression Region	20-122aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	15.8kDa
Protein Length	partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

