



Recombinant Human Multiple inositol polyphosphate phosphatase 1 (MINPP1)


Product Code	CSB-EP891977HU
Relevance	Acts as a phosphoinositide 5- and phosphoinositide 6-phosphatase and regulates cellular levels of inositol pentakisphosphate (InsP5) and inositol hexakisphosphate (InsP6). Also acts as a 2,3-bisphosphoglycerate 3-phosphatase, by mediating the dephosphorylation of 2,3-bisphosphoglycerate (2,3-BPG) to produce phospho-D-glycerate without formation of 3-phosphoglycerate. May play a role in bone development (endochondral ossification).
Abbreviation	Recombinant Human MINPP1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UNW1
Alias	2,3-bisphosphoglycerate 3-phosphatase (EC:3.1.3.80) ;2,3-BPG phosphatase;inositol (1,3,4,5)-tetrakisphosphate 3-phosphatase ;Ins(1,3,4,5)P(4) 3-phosphatase
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	<p> SLLEPRDPVASSLSPYFGTKTRYEDVNPVLLSGPEAPWRDPELLEGTCTPVQL VALIRHGTRYPTVKQIRKLRQLHGLLQARGSRDGGASSTGSRDLGAALADWPL WYADWMDGQLVEKGRQDMRQLALRLASLFPALFSRENYGRLRLITSSKHRC MDSSAAFLQGLWQHYPGLPPPVDADMEFGPPTVNDKLMRFFDHCEKFLTE VEKNATALYHVEAFKTGPEMQNILKKVAATLQVPVNDLNADLIQVAFFTCSDL AIKGVKSPWCDVFDIDDAKVLEYLNDLKQYWKRGYGYTINSRSSCTLFQDIFQ HLDKAVEQKQRSQPISSPVILQFGHAETLLPLLSLMGYFKDKEPLTAYNYKKQM HRKFRSGLIVPYASNLIFVLYHCENAKTPKEQFRVQMLLNEKVLPLAYSQETVS FYEDLKNHYKDILQSCQTSEECELARANSTSEDEL </p>
Research Area	Signal Transduction
Source	E.coli
Target Names	MINPP1
Expression Region	31-487aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged

**Mol. Weight**

79.1kDa

Protein Length

Full Length of Mature Protein

Image

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.