




Recombinant Bordetella pertussis Pertactin autotransporter (prn), partial

Product Code	CSB-EP321224BUA
Relevance	Agglutinin that binds to eukaryotic cells; a process mediated by the R-G-D sequence. Pertactin may have a role in bacterial adhesion, and thus play a role in virulence. May contribute to the disease state of whooping cough.
Abbreviation	Recombinant Bordetella pertussis Pertactin autotransporter protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P14283
Alias	P.93
Product Type	Recombinant Protein
Immunogen Species	Bordetella pertussis (strain Tohama I / ATCC BAA-589 / NCTC 13251)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	ALSKRLGELRLNPDAGGAWGRGFAQRQQLDNRAGRFRFDQKVAGFELGADHA VAVAGGRWHLGGLAGYTRGDRGFTGDGGGHTDSVHVGGYATYIADSGFYLD ATLRASRLENDFKVAGSDGYAVKGYRTHGVGASLEAGRRFTHADGWFLP QAELAVFRAGGGAYRAANGLRVRDEGGSSVLGRLGLEVGKRIELAGGRQVQ PYIKASVLQEFDGAGTVHTNGIAHRTELRGTRAEGLGMAAALGRGHSLYASY EYSKGPKLAMPWTFHAGYRYSW
Research Area	Others
Source	E.coli
Target Names	prn
Expression Region	632-910aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	45.8kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.