




Recombinant Human DnaJ homolog subfamily B member 8 (DNAJB8)

Product Code	CSB-EP836779HU
Relevance	Efficient suppressor of aggregation and toxicity of disease-associated polyglutamine proteins.
Abbreviation	Recombinant Human DNAJB8 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NHS0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MANYEVLGVQASASPEDIKKAYRKLALRWHPDKNPDNKEEAEEKFKLVSEA YEVLSDSKKRSLYDRAGCDSWRAGGGASTPYHSPFDTGYTFRNPEDIFREFF GGLDPFSEFWDSPFNDRGGRGHGLRGAFSAGFGEFPFMEAFSSFNMLG CSGGSHTTFSSTSFSGSSSGSGFKSVMSSTEMINGHKVTTKRIVENGQERV EVEEDGQLKSVTVNGKEQLKWMSK
Research Area	Signal Transduction
Source	E.coli
Target Names	DNAJB8
Expression Region	1-232aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	41.7kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.