



Recombinant Human Heterogeneous nuclear ribonucleoprotein H (HNRNPH1)

Product Code	CSB-EP010609HU
Relevance	This protein is a component of the heterogeneous nuclear ribonucleoprotein (hnRNP) complexes which provide the substrate for the processing events that pre-mRNAs undergo before becoming functional, translatable mRNAs in the cytoplasm. Mediates pre-mRNA alternative splicing regulation. Inhibits, together with CUGBP1, insulin receptor (IR) pre-mRNA exon 11 inclusion in myoblast. Binds to the IR RNA. Binds poly(RG).
Abbreviation	Recombinant Human HNRNPH1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P31943
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MLGTEGGEGFVVKVRGLPWSCSADEVQRFFSDCKIQNGAQQGIRFIYTREGRP SGEAFVELESEDEVKLALKKDRETMGHRVVEVFKSNNVEMDWVLKHTGPNSP DTANDGFVRLRGLPFGCSKEEIVQFFSGLEIVPNGITLPVDFQGRSTGEAFVQF ASQEIAEKALKKKKERIGHRYIEIFKSSRAEVRTHYDPPRKLMMAMQRPQPYDRP GAGRGYNSIGRGAGFERMRRGAYGGGYGGYDDYNGYNDGYGFGSDRFGR DLNYCFSGMSDHRYGDGGSTFQSTTGHCVHMRGLPYRATENDIYNFFSPLNP VRVHIEIGPDGRVTGEADVEFATHEDAVAAMSKDKANMQHRYVELFLNSTAG ASGGAYEHRYVELFLNSTAGASGGAYGSQMMGGMGLSNQSSYGGPASQQL SGGYGGGYGGQSSMSGYDQVLQENS SDFQSNIA
Source	E.coli
Target Names	HNRNPH1
Expression Region	2-449aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	65.1kDa
Protein Length	Full Length of Mature Protein
Image	

