





# Recombinant Mouse Cathepsin S (Ctss)

<b>Product Code</b>	CSB-EP006204MO
<b>Relevance</b>	Thiol protease. Key protease responsible for the roval of the invariant chain from MHC class II molecules. The bond-specificity of this proteinase is in part similar to the specificities of cathepsin L and cathepsin N.
<b>Abbreviation</b>	Recombinant Mouse ctss protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O70370
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	LPDVTVDWREKGCVTEVKYQGSCGACWAFSAVGALEGQLKLTGKLISLSAQN LVDCSNEEKYGNKGC GGGYMTEAFQYIIDNGGIEADASYPYKATDEKCHYNS KNRAATCSRYIQLPFGDEDALKEAVATKGPVSVGIDASHSSFFFYKSGVYDDP SCTGNVNHGVLVVG YGTLDGKDYWLVKNSWGLNFGDQGYIRMARNNKNH GIASYCSYPEI
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	ctss
<b>Expression Region</b>	123-340aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	27.7kDa
<b>Protein Length</b>	Full Length of Mature Protein


## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP006204MO could indicate that this peptide derived from E.coli-expressed Mus musculus (Mouse) ctss.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP006204MO could indicate that this peptide derived from E.coli-expressed Mus musculus (Mouse) ctss.

## Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.