




Recombinant Human NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 6 (NDUFA6)

Product Code	CSB-EP015631HU
Relevance	Accessory subunit of the mitochondrial mbrane respiratory chain NADH dehydrogenase (Complex I), that is believed to be not involved in catalysis. Complex I functions in the transfer of electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone.
Abbreviation	Recombinant Human NDUFA6 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P56556
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MAGSGVRQATSTASTFVKPIFSRDMNEAKRRVRELYRAWYREVPNTVHQFQL DITVKMGRDKVREMFMKNAHVTDPRVVDLLVIKKGIELEETIKVWKQRTTHVMR FFHETEAPRPKDFLSKFYVGHDP
Research Area	Transport
Source	E.coli
Target Names	NDUFA6
Expression Region	1-128 aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	42.1kDa
Protein Length	Full length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.