




Recombinant Human Aquaporin-4 (AQP4), partial

Product Code	CSB-YP001964HU
Relevance	Forms a water-specific channel. Osmoreceptor which regulates body water balance and mediates water flow within the central nervous syst.
Abbreviation	Recombinant Human AQP4 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P55087
Alias	Mercurial-insensitive water channel ;MIWCWCH4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	CPDVEFKRRFKEAFSKAAQQTKGSYMEVEDNRSQVETDDLILKPGVVHVIDVD RGEEKKGKDKSGEVLSSV
Research Area	Transport
Source	Yeast
Target Names	AQP4
Protein Names	Recommended name: Aquaporin-4 Short name= AQP-4Alternative name(s): Mercurial-insensitive water channel Short name= MIWC WCH4
Expression Region	253-323aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	10.0kDa
Protein Length	Cytoplasmic Domain

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.