




# Recombinant Human DENN domain-containing protein 1A (DENND1A)

<b>Product Code</b>	CSB-EP823472HU
<b>Relevance</b>	Guanine nucleotide exchange factor (GEF) regulating clathrin-mediated endocytosis through RAB35 activation. Promotes the exchange of GDP to GTP, converting inactive GDP-bound RAB35 into its active GTP-bound form. Regulates clathrin-mediated endocytosis of synaptic vesicles.
<b>Abbreviation</b>	Recombinant Human DENND1A protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8TEH3
<b>Alias</b>	Connecdenn 1 ;Connecdenn;Protein FAM31A
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 90% as determined by SDS-PAGE.
<b>Sequence</b>	MGSRIKQNPETTFEVYVEVAYPRTGGTSLSDPEVQRQFPEDYSDQEVLTQTLTKF CFPFYVDSLTVSQVGQNFVLTLDIDSKQRFGFCRLSSGAKSCFCILSYLPWFE VFYKLLNILADYTTKRQENQWNELLETLHKLPIPDGVSVHLSVHSYFTVPDTR ELPSIPENRNLTEYFVAVDVNNMLHLYASMLYERRILIICSKLSTLTACIHGSAAM LYPMYWQHVIYIPVLPPLLDYCCAPMPYLIIGIHLSLMEKVRNMLDDVVILNVD TNTLETPFDDLQSLPNDVISSLKNRLKVVSTTTGDGVARAFLKAQAFFGYSYRN ALKIEPEEPIFCEEAFVSHYRSGAMRQFLQATQLQLFKQFIDGRLDLLNSGE GFSDVFEIEINMGEYAGSDKLYHQWLSTVRKGS GAILNTVKTKANPAMKTVYK FAKDHA KMGIKEVKNRLKQKDIAENG CAPTPEEQLPKTAPSPLVEAKDPKLRE DRRPITVHFGQVRPPRPHVVKRPKSNIAVEGRRTSVPSPEQNTIATPATLHILQ KSITHFAAKFPTRGWTSSSH
<b>Research Area</b>	Neuroscience
<b>Source</b>	E.coli
<b>Target Names</b>	DENND1A
<b>Expression Region</b>	1-559aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-SUMO-tagged
<b>Mol. Weight</b>	79.4kDa
<b>Protein Length</b>	Full Length of Isoform 2



## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

## Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.