



Recombinant Human Guanine nucleotide-binding protein-like 3-like protein (GNL3L)


Product Code	CSB-EP865150HU
Relevance	Stabilizes TERF1 telomeric association by preventing TERF1 recruitment by PML. Stabilizes TERF1 protein by preventing its ubiquitination and hence proteasomal degradation. Does so by interfering with TERF1-binding to FBXO4 E3 ubiquitin-protein ligase. Required for cell proliferation. By stabilizing TRF1 protein during mitosis, promotes metaphase-to-anaphase transition. Stabilizes MDM2 protein by preventing its ubiquitination, and hence proteasomal degradation. By acting on MDM2, may affect TP53 activity. Required for normal processing of ribosomal pre-rRNA. Binds GTP.
Abbreviation	Recombinant Human GNL3L protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9NVN8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MMKLRHKNKKPGEGSKGHKKISWPYPQPAKQNGKKATSKVPSAPHFVHPND HANREAELKKKWVEEMREKQQAAREQERQKRRTIESYCQDVLRRQEEFEHK EEVLQELNMFQQLDDEATRKAYYKEFRKVVEYSDVILEVLDARDPLGCRCFQM EEAVLRAQGNNKKLVVLNKKIDLVPEVVEKWL DYLRNELPTVAFKASTQHVK NLNRCVPVDQASESLKSKACFGAENLMRVLGNYCRLGEVRTHIRVGVVGL PNVGKSSLINSLKRSRACSVGAVPGITKFMQEVYLDKFIRLLDAPGIVPGPNSE VGTILRNCVHVQKLADPVTPTVETILQRCNLEEISNYYGVSGFQTTEHFLTAVAH RLGKKKKGGLYSQEQAAKAVLADWVSGKISFYIPPPATHTLPTHL SAEIVKEMT EVFDIEDTEQANEDTMECLATGESDELLGDTDPLEMEIKLLHSPMTKIADAIENK TTVYKIGDLTGYCTNPNRHQMGWAKRNVDHRPKSNSMVDVCSDRRSVLQR IMETDPLQQGQALASALKNKKKMQRADKIASKLSDSMMSALDLSGNADDGV GD
Research Area	Others
Source	E.coli
Target Names	GNL3L
Expression Region	1-582aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged



Mol. Weight 81.6kDa

Protein Length Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.