



MRPL42 Antibody

Product Code	CSB-PA014852EA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6G3
Immunogen	Recombinant Human 39S ribosomal protein L42, mitochondrial protein (33-142AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	39S ribosomal protein L42, mitochondrial (L42mt) (MRP-L42) (39S ribosomal protein L31, mitochondrial) (L31mt) (MRP-L31) (Mitochondrial large ribosomal subunit protein mL42), MRPL42, MRPL31 MRPS32 RPML31
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	MRPL42
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.