



Recombinant Human HLA class I histocompatibility antigen, alpha chain G (HLA-G)

Product Code	CSB-MP010509HU
Relevance	Involved in the presentation of foreign antigens to the immune syst. Plays a role in maternal tolerance of the fetus by mediating protection from the deleterious effects of natural killer cells, cytotoxic T-lymphocytes, macrophages and mononuclear cells.
Abbreviation	Recombinant Human HLA-G protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P17693
Alias	HLA G antigen;MHC class I antigen G
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	GSHSMRYFSAAVSRPGRGEPRIAMGYVDDTQFVRFDSDSACPRMEPRAPW VEQEGPEYWEEETRNTKAHAQTDRMNLQTLRGYYNQSEASSHTLQWMIGCD LGSDGRLLRGYEQYAYDGKDYALNEDLRSWTAADTAAQISKRKCEANVAE QRRAYLEGTCVEWLHRYLENGKEMLRADPPKTHVTHHPVFDYEATLRCWA LGFYPAEIIITWQRDGEDQTQDVELVETRPAGDGTGFQKWAAVVVPSGEEQRY TCHVQHEGLPEPLMRWKQSSLPTIPIMGIVAGLVVLAAVVTGAAVA AVLWRK KSSD
Source	Mammalian cell
Target Names	HLA-G
Expression Region	25-338aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	39.6kDa
Protein Length	Full Length of Mature Protein
Image	

