



A1BG Antibody

Product Code	CSB-PA001001GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04217
Immunogen	Human A1BG
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	alpha-1-B glycoprotein;A1BG;A1B;ABG;DKFZp686F0970;GAB;HYST2477 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	A1BG
Target Details	This protein is a plasma glycoprotein of unknown function. The protein shows sequence similarity to the variable regions of some immunoglobulin supergene family member proteins.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.