



AAAS Antibody

Product Code	CSB-PA001010GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRG9
Immunogen	Human AAAS
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	achalasia, adrenocortical insufficiency, alacrimia (Allgrove, triple-A);AAAS;AAA;AAASb;ADRACALA;ADRACALIN;ALADIN;DKFZp586G1624;GLO03 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	AAAS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.