



ABHD14A Antibody

Product Code	CSB-PA001091GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BUJ0
Immunogen	Human ABHD14A
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	abhydrolase domain containing 14A;ABHD14A;DKFZp564O243;DORZ1;FLJ60317;FLJ60467 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ABHD14A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.