



ABT1 Antibody

Product Code	CSB-PA001114GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULW3
Immunogen	Human ABT1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	activator of basal transcription 1;ABT1;hABT1 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ABT1
Target Details	Basal transcription of genes by RNA polymerase II requires the interaction of TATA-binding protein (TBP) with the core region of class II promoters. Studies in mouse suggest that This protein likely activates basal transcription from class II promoters by interaction with TBP and the class II promoter DNA.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.