



ACSS1 Antibody

Product Code	CSB-PA001202GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NUB1
Immunogen	Human ACSS1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	acyl-CoA synthetase short-chain family member 1;ACSS1;ACAS2L;AceCS2L;FLJ45659;MGC33843 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ACSS1
Target Details	An acetyl-CoA synthetase found in the mitochondrial matrix has been characterized in both cow and mouse. It is primarily a cardiac enzyme which produces acetyl-CoA mainly for the oxidation of acetate. This protein is structurally similar to the cow and mouse proteins and therefore may serve a similar function in humans.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.