



ACYP2 Antibody

Product Code	CSB-PA001266GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14621
Immunogen	Human ACYP2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	acylphosphatase 2, muscle type;ACYP2;ACYM;ACYP ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ACYP2
Target Details	Acylphosphatase can hydrolyze the phosphoenzyme intermediate of different membrane pumps, particularly the Ca ²⁺ /Mg ²⁺ -ATPase from sarcoplasmic reticulum of skeletal muscle. Two isoenzymes have been isolated, called muscle acylphosphatase and erythrocyte acylphosphatase on the basis of their tissue localization. This gene encodes the muscle-type isoform (MT). An increase of the MT isoform is associated with muscle differentiation.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.