



AGPAT9 Antibody

Product Code	CSB-PA001456GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q53EU6
Immunogen	Human AGPAT9
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IF
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	1-acylglycerol-3-phosphate O-acyltransferase 9;AGPAT9;AGPAT8;GPAT3;HMFN0839;LPAAT-theta;MAG1;MGC11324 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	GPAT3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.