



PHYKPL Antibody

Product Code	CSB-PA001471GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IUZ5
Immunogen	Human AGXT2L2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	alanine-glyoxylate aminotransferase 2-like 2;AGXT2L2;MGC117348;MGC15875;MGC45484 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PHYKPL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.