



ALKBH7 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA001613GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9BT30 |
| Immunogen | Human ALKBH7 |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,WB |
| Storage Buffer | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity purified |
| Isotype | IgG |
| Alias | alkB, alkylation repair homolog 7 (E. coli);ALKBH7;MGC10974;SPATA11;UNQ6002 ; |
| Product Type | Purified Rabbit Anti human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | ALKBH7 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |