



AMICA1 Antibody

Product Code	CSB-PA001668GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86YT9
Immunogen	Human AMICA1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	adhesion molecule, interacts with CXADR antigen 1;AMICA1;AMICA;CREA7-1;CREA7-4;FLJ37080;Gm638;JAML;MGC118814;MGC118815 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	JAML
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.