



ANP32A Antibody

Product Code	CSB-PA001822GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P39687
Immunogen	Human ANP32A
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	acidic (leucine-rich) nuclear phosphoprotein 32 family, member A;ANP32A;C15orf1;I1PP2A;LANP;MAPM;MGC119787;MGC150373;PHAP1;PHAPI;PP32 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ANP32A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.