



ANTXR2 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA001833GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P58335 |
| Immunogen | Human ANTXR2 |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | ELISA,IHC |
| Storage Buffer | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity purified |
| Isotype | IgG |
| Alias | anthrax toxin receptor 2;ANTXR2;CMG-2;CMG2;FLJ31074;ISH;JHF;MGC111533;MGC45856 ; |
| Product Type | Purified Rabbit Anti human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | ANTXR2 |
| Target Details | This gene encodes a receptor for anthrax toxin. The protein binds to collagen IV and laminin, suggesting that it may be involved in extracellular matrix adhesion. Mutations in this gene cause juvenile hyaline fibromatosis and infantile systemic hyalinoses. Multiple transcript variants encoding different isoforms have been found for this gene. |