



APLF Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA001908GA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8IW19 |
| Immunogen | Human APLF |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, WB |
| Storage Buffer | PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles. |
| Purification Method | Antigen Affinity purified |
| Isotype | IgG |
| Alias | aprataxin and PNKP like factor; APLF; APFL; C2orf13; FLJ16593; MGC47799; PALF; Xip1 ; |
| Product Type | Purified Rabbit Anti human PolyClonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | APLF |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |