



ARHGAP24 Antibody

Product Code	CSB-PA002025GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N264
Immunogen	Human ARHGAP24
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat, Dog, Zebrafish
Tested Applications	ELISA, WB, IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	Rho GTPase activating protein 24; ARHGAP24; DKFZp564B1162; FILGAP; FLJ33877; RC- GAP72; RCGAP72; p73; p73RhoGAP ; FILGAP
Product Type	Purified Rabbit Anti Human, Mouse PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ARHGAP24
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.