



ARHGDIB Antibody

Product Code	CSB-PA002038GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52566
Immunogen	Human ARHGDIB
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	Rho GDP dissociation inhibitor (GDI) beta;ARHGDIB;D4;GDIA2;GDID4;LYGDI;Ly-GDI;RAP1GN1;RhoGDI2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ARHGDIB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.