



ARL6IP1 Antibody

Product Code	CSB-PA002094GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15041
Immunogen	Human ARL6IP1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	ADP-ribosylation factor-like 6 interacting protein 1;ARL6IP1;AIP1;ARL6IP;ARMER;KIAA0069 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ARL6IP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.