



ASH2L Antibody

Product Code	CSB-PA002210GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBL3
Immunogen	Human ASH2L
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	ash2 (absent, small, or homeotic)-like (Drosophila);ASH2L;ASH2;ASH2L1;ASH2L2;Bre2 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ASH2L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.