



ASNA1 Antibody

Product Code	CSB-PA002218GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43681
Immunogen	Human ASNA1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	arsA arsenite transporter, ATP-binding, homolog 1 (bacterial);ASNA1;ARSA-I;ARSA1;MGC3821 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ASNA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.