



ATG16L1 Antibody

Product Code	CSB-PA002284GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q676U5
Immunogen	Human ATG16L1
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity purified
Isotype	IgG
Alias	ATG16 autophagy related 16-like 1 (S. cerevisiae);ATG16L1;APG16L;ATG16L;FLJ00045;FLJ10035;FLJ10828;FLJ22677;IBD10;WDR30 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ATG16L1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.