



ATP6V1E2 Antibody

Product Code	CSB-PA002403GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96A05
Immunogen	Human ATP6V1E2
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	ATPase, H+ transporting, lysosomal 31kDa, V1 subunit E2;ATP6V1E2;ATP6E1;ATP6EL2;ATP6V1EL2;MGC9341;VMA4 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ATP6V1E2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.