



BAP1 Antibody

Product Code	CSB-PA002556GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92560
Immunogen	Human BAP1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB
Storage Buffer	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Alias	BRCA1 associated protein-1 (ubiquitin carboxy-terminal hydrolase); BAP1; DKFZp686N04275; FLJ35406; FLJ37180; HUCEP-13; KIAA0272; UCHL2; hucep-6 ;
Product Type	Purified Rabbit Anti human PolyClonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	BAP1
Target Details	This protein localizes to the nucleus and it interacts with the RING finger domain of the breast cancer 1, early onset protein (BRCA1). This gene is thought to be a tumor suppressor gene that functions in the BRCA1 growth control pathway. There are multiple polyadenylation sites found in this gene.
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.